

Patient information leaflet for KHEIRON Spinal Fixation System

1. What is KHEIRON Spinal Fixation System used for

The KHEIRON Spinal Fixation System including the K-ROD patient-specific rod is a spinal fixation system intended to help provide immobilisation and stabilisation of spinal segments as an adjunct to bone fusion of the thoracic, lumbar, and/or sacral spine.

These devices are indicated in skeletally mature patients when conservative treatments were not efficient or when the disease progression may represent a threat to patient safety, for all the following indications:

- Degenerative disc disease (defined as back pain with degeneration of the disc confirmed by history and radiographic studies),
- Spondylolisthesis, i.e sliding of a vertebra forward and downward in relation to the vertebrae just below it,
- Trauma (i.e., fracture or dislocation),
- Deformities or curvatures (i.e., scoliosis, kyphosis, and/or lordosis, Scheuermann's Disease),
- Tumour,
- Stenosis (i.e narrowing of the spinal canal)
- Pseudoarthrosis (i.e failed previous bone fusion).
- Treatment of severe spondylolisthesis of the L5-S1 vertebrae.

2. Product description

KHEIRON Spinal Fixation System includes pedicle screws, straight and curved rods of Ø5.5mm and Ø6mm, transverse systems, dominos connectors, sacral plates, sacral screws, iliac extension, iliac screws, hooks, claw hooks and connectors. Those components can be locked in various configurations, each assembly being tailor-made.

K-ROD rods are manufactured for a specific patient, each rod designed is identified by a patient code (ID KEOPS). The K-ROD patient specific rods are available in Ø5.5 and Ø6mm.

All implants are made of titanium alloy.

3. When KHEIRON Spinal Fixation System should not be used

Contraindications include, but are not limited to:

- The bony condition of the patient (e.g.: massive or severe osteoporosis) making the procedure risky in terms of mechanical securing of the implant
- Congenital spinal stenosis
- Comminuted (multiple) fractures involving several vertebrae
- Tumours involving several successive vertebrae
- Allergies, intolerance and/or hypersensitivity to the component material of the implant which is Ti-6Al-4V ELI, titanium alloys
- Primary or secondary infection
- Local inflammation
- Fever, leukocytosis
- Obesity
- Pregnancy
- Mental illness or patient not likely to follow the surgeon's recommendations
- Congenital abnormal anatomy
- Rapidly progressive joint disease, severe osteoporosis
- Inappropriate anatomy
- Patient having inadequate tissue coverage over the operative site

Patient information leaflet for KHEIRON Spinal Fixation System

- Any case not described in the indications described in section 1.

4. Caution

The following information should be considered before the procedure in order to ensure success of the surgical implantation:

- Clinical data show that patients who smoke tend to have less optimum bony consolidation, as do patients who are undernourished or alcoholic.
- To aid bone healing, it is important to limit use of nicotine and non-steroidal medicinal products (e.g.: aspirin).
- The implanted device must not be exposed to unwanted forces such as mechanical vibrations. Consequently, the patient must limit his or her physical activity (athletic and occupational), especially in cases of lifting, twisting and crushing.
- Throughout the consolidation period, the patient must follow the surgeon's instructions and recommendations.

5. Adverse effects

The following adverse effects have been observed. This list may not include all complications caused by the surgical technique itself:

- Displacement or expulsion of the implant before bone fusion requiring another surgical procedure
- Device malposition potentially requiring revision surgery
- Haematoma
- Infection at implantation site
- Failed bone fusion (pseudarthrosis)
- Neurological damage
- Bone damage
- Disassembly, deformation and/or rupture of one or more components of the system
- Risk of allergy to titanium alloy is rare, but must be considered
- Interference with radiographic, CT scan, and/or MR imaging caused by the implants
- Complications due to surgery (haemorrhage, infection, complications due to the use of bone grafting, respiratory problems, reaction to anaesthesia, death).

6. MRI information

- These implants do not present any known risk of interference with other medical equipment.
- Device safety and compatibility in a magnetic resonance setting have not been evaluated. No thermal test or migration test has been performed on the device in this setting.

7. Symptoms that could be related to a dysfunction of the device

Fever, back pain and squeaks at the operated site may be the signs of a dysfunction of the device.

8. Lifetime

The KHEIRON Spinal Fixation System is a spinal fixation system intended to help provide immobilisation and stabilisation of spinal segments as an adjunct to bone fusion. Bone fusion may occur after a variable period depending on each patient. After fusion occurs, the material can be removed (in this case the benefit/risk ratio of an additional surgery should be discussed with the practitioner).

Patient information leaflet for KHEIRON Spinal Fixation System

9. Monitoring/maintenance of the implanted devices

No specific monitoring/maintenance of the implanted devices is necessary. In case of adverse event your surgeon may request a radiographic check-up.

10. Vigilance

Any adverse outcomes potentially attributable to KHEIRON Spinal Fixation System implants must be reported promptly to your doctor. You may also report it to S.M.A.I.O at vigilance@smaio.com or directly to the Therapeutic Goods Administration at this link: <http://www.tga.gov.au/reporting-problems>

11. Manufacturer

S.M.A.I.O
 2, Place Berthe Morisot – Parc Technologique
 69800 Saint-Priest
 FRANCE
 Tel.: +33 (0)4 69 84 23 02
 Website: www.smaio.com

12. List of devices covered by this leaflet

24BR55XXX	26PH00001	27AC51000	28IE56XXX	29TC55XXX	32KR55000
24BR60XXX	26PH45XXX	27PC52XXX	28IS00001	29TC60XXX	32KR60000
24CR55XXX	26PH67XXX	27PC53XXX	28IWXXXXXX		
24CR60XXX	26SL00XXX	27PC54XXX			
24PW55600	26SP45XXX	27PC55XXX			
24PW55700	26SP67XXX				
24PW56XXX	26TC00001				
24PW60600	26TC10XXX				
24PWXXXXX	26TC20XXX				
24SR55XXX					
24SR60XXX					
25SC55XXX					
25SC60XXX					
25SN55XXX					
25SP55XXX					
25SWXXXXX					

Where “XXX” and “XXXXX” can be any combination of numbers.